

Technical Data Sheet



1. Product


050 - pH Indicator Paper, 0.5 - 5.0

2. Instructions for use

- Immerse the indicator paper in the test solution for 1 second.
- Compare the colour produced on the strip against the colour chart before the indicator paper dries out.

3. Measuring range

pH 0.5 - 1.0 - 1.5 - 2.0 - 2.5 - 3.0 - 3.5 - 4.0 - 4.5 - 5.0



4. Application

pH indicator paper for testing of pH in solutions between pH 0.5 and pH 5.5.

5. Shelf Life

3 years from manufacturing date.

6. Contents

100 strips.

7. Weights & Dimensions

pH paper: 5m x 7mm

8. Notes

Keep any unused strips in the original container. Do not store where they are liable to be exposed to acid or alkaline gases or vapours. Store in a cool dry place.

Johnson Test Papers Ltd
Hainge Park, Hainge Road
Oldbury, West Midlands
B69 2NU, United Kingdom

T: +44 (0)121 557 3883
F: +44 (0)121 557 8235
E: info@jtp.uk.com
W: www.johnsonstestpapers.com

Registered in England
Company Reg No: 4112470
VAT No: GB 754 6913 04

