

# Technical Data Sheet

## 1. Product

224.1 - J-QUANT Phosphate 0 - 500 mg/L, 100 strips/pk

## 2. Instructions for use

- Immerse an indicator strip in the test solution for 10 seconds.
- Shake off excess solution.
- Wait 3 minutes.
- Compare the colour produced on the strip against the colour chart.

## 3. Measuring range

0 - 3 - 10 - 25 - 50 - 100 - 150 - 500 mg/L



## 4. Application

The J-QUANT Phosphate test strip detects the level of orthophosphate ions in solution, for example, in waste water, soil samples, fertilizers and food. They can also be used for measuring corrosion inhibitors when servicing cooling and heating systems.

## 5. Shelf Life

3 years from manufacturing date.

## 6. Contents

100 strips.

## 7. Weights & Dimensions

Strip dimension: 83mm x 5mm.

## 8. Notes

Keep any unused strips in the original container. Do not store where they are liable to be exposed to acid or alkaline gases or vapours. Store in a cool dry place.