

Technical Data Sheet



1. Product

172.5C - J-QUANT Sulphite 0 - 500 mg/L, 100 strips/pk

2. Instructions for use

- Ensure the pH of the solution to be tested is between 8 and 9. Adjust if necessary.
- Immerse an indicator strip in the test solution for 1 second.
- Shake off excess solution.
- Wait 10 seconds.
- Compare the colour produced on the strip against the colour chart.

3. Measuring range

0 - 10 - 25 - 50 - 100 - 250 - 500 mg/L



4. Application

The J-QUANT Sulphite test strip detects the level of sulphite, for example, in boiler water, fruit juices and the wine making process. This test strip is compatible with all forms of sulphite ion (sulphite, bisulphite and metabisulphite).

5. Shelf Life

2 years from manufacturing date.

6. Contents

100 strips

7. Weights & Dimensions

Strip dimension: 83mm x 5mm

8. Notes

Keep any unused strips in the original container. Do not store where they are liable to be exposed to acid or alkaline gases or vapours. Store in a cool dry place.

Johnson Test Papers Ltd
Hainge Park, Hainge Road
Oldbury, West Midlands
B69 2NU, United Kingdom

T: +44 (0)121 557 3883
F: +44 (0)121 557 8235
E: info@jtp.uk.com
W: www.johnsonstestpapers.com

Registered in England
Company Reg No: 4112470
VAT No: GB 754 6913 04



ISO 9001 : 2015



0273



ISO 13485 : 2014



0273