

# Technical Data Sheet

## 1. Product

218.1 - J-QUANT Nitrite 0 - 80 mg/L, 100 strips/pk

## 2. Instructions for use

- Immerse an indicator strip in the test solution for 1 second.
- Shake off excess solution.
- Wait 45 seconds.
- Compare the colour produced on the pH strip with the colour chart.

## 3. Measuring range

0 - 2 - 5 - 10 - 20 - 40 - 80 mg/L



## 4. Application

The J-QUANT Nitrite test strip detects residual and build-up of nitrite, for example, in drinking water, waste water or aquariums.

## 5. Shelf Life

2 years from manufacturing date.

## 6. Contents

100 strips.

## 7. Weights & Dimensions

Strip dimension: 83mm x 5mm.

## 8. Notes

Keep any unused strips in the original container. Do not store where they are liable to be exposed to acid or alkaline gases or vapours. Store in a cool dry place.